

Introducing... **Dragonfly**

The most flexible imaging solution ever!

1-2 marzo

The Complete Solution.

The **Dragonfly** naturally benefits from Andor EMCCD and sCMOS technologies, with motorized optical zoom, the system delivers industry-leading signal to noise and image fidelity.

Applications ranging from single molecule to live cell confocal, TIRFM to whole embryo and thick tissue imaging all benefit from the system's speed and sensitivity.

Dragonfly delivers real-time visualization for rapid specimen evaluation and GPU-accelerated deconvolution to maximize resolution and throughput.

- Instant Confocal
- Simultaneous dual color TIRF
- Laser widefield imaging
- Borealis Illumination
- NIR extended range

UNITAT DE
MICROSCOPIA ÒPTICA AVANÇADA
CCiT UB CAMPUS CLÍNIC
FACULTAT DE MEDICINA
C/ CASANOVAS, 143
ALA SUD 5º PLANTA
BARCELONA



 **Izasa
Scientific**

A WerfenLife Company

izasascientific.com

ANDOR
an Oxford Instruments company